Cooper Surgical

Wallach® WA1000B™ General Purpose Cryosurgical System (LL-CO₂™)
900162, 900163, 900509-5, 900509-84, 900509-184,
900506-1-CO₂, 900506-2-CO₂

Instructions for Use (English)

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CE₂₇₉₇

The CE mark does not apply to REF Numbers 900509-5, 900509-184, 900506-1-CO₂ and 900506-2-CO₂.

Instructions for Use (English)	Page 1
Gebruiksaanwijzing (Nederlands / Dutch)	Pagina 8
Instrucciones de uso (Español [América Latina] / Spanish [Latin America]).	Página 15
Petunjuk Penggunaan (Bahasa Indonesia / Indonesian)	Halaman 22
คำแนะนำการใช้งาน (ภาษาไทย / Thai)	หน้า 29

**Wallach® WA1000B™ General Purpose Cryosurgical System (LL-CO₂™)
900162, 900163, 900509-5, 900509-84, 900509-184,
900506-1-CO₂, 900506-2-CO₂**

Instructions for Use (English)

READ ALL SAFETY INFORMATION AND INSTRUCTIONS BEFORE USING THIS PRODUCT!

DEVICE DESCRIPTION

The Wallach® WA1000B™ Cryosurgical System (LL-CO₂) is a product of CooperSurgical. It consists of a Cryo Console and Cryosurgical Tips. The system requires the use of Medical Grade CO₂ (Carbon Dioxide) gas.

WA1000B Cryo Console

The WA1000B Cryo Console consists of a freezer handle and probe with attached tubing. It is also provided with a stand and a cylinder yoke connector.

Cryosurgical Tips

A variety of reusable cryotips with disposable plastic shields are available for use with the unit.

INDICATIONS FOR USE

The Wallach WA1000B Cryosurgical System (LL-CO₂) is used for a variety of gynecology and dermatology cryosurgical procedures in the process of controlled destruction of defined areas of benign or malignant lesions.

WARNINGS

- **Do not use overfilled cylinder; gas pressure may be excessive and present a hazard. If gas pressure gauge reads below 45.70 kg/cm², replace the cylinder. If cylinder pressure is above 56.25 kg/cm², bleed excess gas from cylinder. Do not operate freezer if gauge is in the red zone.**
- **Visually inspect the cryo instrument and tips prior to each use. Discontinue use of any product suspected of being damaged. Such product should be returned to CooperSurgical for inspection and service.**
- **At no time should chlorine mixtures (bleach) or iodine-based products be used for soaking or cleaning Wallach cryotips. Cryotips will degrade if exposed to such solution.**
- **In case of malfunction, close gas cylinder master valve.**
- **In the event of a defrost failure, turn tank off; leave tip in place. Remove tip when defrost is evident or use warm water to expedite.**
- **Changing tips must be done with console switch “OFF”.**
- **Do not soak cryotips as lines can become clogged.**
- **Moisture in lines and orifices of the cryotips and cryo console can cause sputtering, poor freeze and other malfunctions. Whenever using fluids or steam to clean or sterilize, all gas lines and cryotip orifices must be plugged.**

CAUTION

- **U.S. Federal law restricts this device to sale by or on the order of a physician.**
- **At the end of each procedure or when changing tips make sure console switch is turned off. Freezer trigger should be flipped to the upward position so gas will release through the exhaust tube. Removal of the cryotip prior to 5 minutes after shutdown may result in unexpected cryo damage to unprotected epithelium of the operator’s hand.**
- **Operating cylinder must be at room temperature in order to stabilize the gas pressure (68 °F to 74 °F, 18 °C to 22 °C). Cold will decrease the system’s efficiency and significantly increase defrost time. Heat causes excess gas consumption.**
- **Allow time for tank temperature to recover to room temperature. Do not store gas cylinders near heaters, radiators, or other heat sources, as heat will increase the gas pressure within the cylinder.**
- **Move cylinder to a well-ventilated area when venting CO₂. The gas exhaust port will accept an easily attached scavenger hose.**
- **Do not mix or use gases other than CO₂.**
- **Do not use wrenches, hammers, or other tools to connect Wallach equipment to gas cylinders.**
- **Use only non-siphon type gas cylinders.**
- **To ensure optimal performance, only medical grade gas free of contaminants should be used.**

WA1000B™ CRYOSURGICAL SYSTEM (LL-CO₂) COMPONENTS

1. GAS CYLINDER PRESSURE Gauge	Once cylinder is turned “ ON ” this gauge should read in the 45.70 to 56.25 kg/cm ² range.
2. CONSOLES a) CONSOLE SWITCH b) CONSOLE RECEPTACLE	a. When turned to the “ ON ” position permits gas flow into the probe. b. Freezer connects to the console by screwing the knurled freezer connector to the receptacle.
3. GAS EXHAUST PORT	A gas scavenger hose may be attached to this port.
4. THE TIMER/TEMPERATURE INDICATOR CONSOLE (WA2000)	Monitors: • Time (in seconds) during which the probe tip experiences freezing temperatures. • Temperature of the probe tip, when used with probes which have an integral thermocouple. The Console will accept an integral thermocouple. You will be able to obtain a digital temperature reading of tip temperature. NOTE: 4 “D” cell batteries are needed for this console. (Console will display “Low Batteries” when batteries are in need of replacement.)
5. LL-CO ₂ ™ Freezer	Multi-tip freezer plugs into console.
6. POWER Switch	Used to turn console power “ ON ” or “ OFF ”.
7. SELECTOR Switch	Used for determining which of two thermocouple ports will be shown on the visual temperature readout panel.
8. TEMPERATURE INDICATOR Display	Readout panel for tip temperature.

9. TIMER Display	The timer segment monitors the time (in seconds) during which the probe tip experiences temperatures of -20°C or colder. Once the probe tip reaches -20°C , the timer will be activated and it will emit audible beeps at 15 second intervals, thus aiding the practitioner in determining the total freezing time. As the probe tip begins to freeze, the timer automatically activates the digital counter, and the timer continues to run until the probe begins a defrost mode, warming up to 0°C , at which point the timing stops and the digital readout holds its current value. The practitioner has the option of resetting the timer to “zero” for a new count, or continuing the current count (cumulative) with an additional freeze. To reset to “zero” at any time, depress the “Reset” switch once.
10. Mobile Cart (WA1500)	Compact cabinet secures and conceals gas cylinder and allows transport of unit.

ASSEMBLY INSTRUCTIONS

1. Attach the two console units to the cart by aligning the holes in the bottom of each unit with the locating pins on the top of the cart. Press each console down to firmly seat it in place on the cart. (The Timer/Temperature Indicator Console mounts on the upper rear part of the cart and the General Purpose Probe Console mounts on the upper front part of the cart.) To remove the console from the cart, simply pull straight up and the consoles will detach themselves from the cart.
2. Insert the stainless steel freezer holder into the hole indicated on the General Purpose Probe Console.
3. To install a gas cylinder in the cart, separate the Velcro straps at the rear of the cart, place the cylinder in position (with the hook-up valve facing out) firmly against the cart frame and tightly secure the straps. Attach the gas connector hose of the Probe Console to the gas cylinder. Hand tighten only.
NOTE: Use only “medical grade” Carbon Dioxide in Non-Siphon Gas cylinders.
4. Connect the LL-CO₂™ Freezer to the CONSOLE RECEPTACLE. Push the connector in and secure it in place by tightening the knurled ring of the probe connector. (Make it snug; hand tighten only.) If a freezer with an integral thermocouple is selected, attach the thermocouple plug to the thermocouple receptacle (TC1 or TC2) on the Timer/Temperature Indicator Console. Be sure to match the color code of the thermocouple prongs to that of the receptacle, in order to assure a proper fit. Move the selector switch to the correct position for the thermocouple receptacle selected.
5. With the console switch in the **OFF** “O” position, select the cryotip to be used. The cryotip need only be finger tight. Do not use tools. Be sure cryotip shield is placed on the tip so that the grooves fit onto the notches on the cryotip. This will assure not only proper seating of the shield on the cryotip but also will provide a positive stop to ensure the cryotip is firmly secured to the probe.
6. To reduce the concentration of CO₂ in the immediate cryosurgery area, you may vent the exhaust gas by attaching one end of the scavenger hose to the GAS EXHAUST PORT. The other end of the hose goes to a discharge vent.

INSTRUCTIONS FOR USE

1. Slowly open the gas cylinder master valve by turning the hand wheel counter-clockwise. Listen for any gas leaks. (There should be none.) Proceed to open the cylinder valve completely.
2. Once the gas pressure is stable, turn CONSOLE SWITCH to the “**ON**” position. The Freezer is now ready to operate.

PRECAUTION: Make sure locking/freeze trigger is not engaged before turning console switch to the “**ON**” position.

NOTE: It is recommended that an agent be used to enhance lubrication and moisture between the freezer and tissue. The surface of the probe should be positioned so as to face the surface of the tissue.

3. Depress the freeze trigger of the Wallach® LL-CO₂™ Freezer. Within seconds you will observe frost on the cryotip. To defrost, release the trigger and the frost will disappear. The freeze trigger will lock in position for a longer freeze if pushed forward. To defrost release trigger. You will see the frost disappear almost immediately. The sound of escaping gas is the normal gas exhaust coming out of the GAS-EXHAUST PORT. (As mentioned, this gas may be vented via an optional scavenger hose.)

NOTE: After a prolonged freeze, the white silicone hose exiting from the freezer handle may become cold. A light frosting may occur, but this is normal. Though cold to the touch, the hose will not freeze or damage tissue.

4. Within 3 seconds after releasing the trigger, the instrument is ready to refreeze.
5. When the procedure has been completed, shut the gas cylinder master valve “**OFF**” by turning the hand wheel clockwise.
6. Turn the CONSOLE SWITCH “**OFF**”. Push the Freeze trigger in the forward position so excess gas will release through exhaust port.

CAUTION: Unit must be completely shut down to change tips during a procedure. Wait 5 minutes before changing the tip to prevent unexpected cryo damage to unprotected epithelium.

7. You are now ready to remove the LL-CO₂ Freezer. If a freezer with an integral thermocouple has been used, remove the probe thermocouple plug from the thermocouple receptacle. Remove the LL-CO₂ Freezer by loosening the knurled ring on the probe connector and pulling the probe connector gently, straight out of the CONSOLE RECEPTACLE. **DO NOT PULL ON THE FLEXIBLE HOSE.** Always use the body of the LL-CO₂ Freezer connector when connecting or disconnecting the probe.
8. Remove the gas connector hose of the Probe Console from the gas cylinder.

TROUBLESHOOTING

PROBLEM	POSSIBLE CAUSE	CORRECTIVE ACTION
1. Gas Cylinder Pressure gauge reads in the RED	The cylinder is over pressurized.	Take immediate action to reduce the gas pressure by closing the cylinder valve, loosening the Cryo connector to the cylinder valve, and carefully turning the cylinder valve “ ON ” again to bleed off the excess pressure. This procedure is best performed in a wide-open or outdoor area (in order to disperse the gas).
2. Gas Cylinder Pressure gauge reads in the YELLOW	There is not enough gas pressure to start a procedure.	Either change the cylinders, or, if the cylinder was stored in a cold area, allow the cylinder to warm up to ambient room temperature
3. Insufficient Freeze to Probe Tip	<ol style="list-style-type: none"> a. Gas leaking from device b. Non-Medical grade gas used c. Functional failure d. Repair done by unqualified personnel e. Operating in the yellow zone on the pressure gauge 	<ol style="list-style-type: none"> a-d. Refer to the service and repair section located in the instruction manual. e. Replace tank
4. Pressure Gauge Failure	Pressure lines clogged inside console.	Purge the gas lines to make sure that they are clear of moisture or foreign matter.
5. No Timer/Temperature Display	<ol style="list-style-type: none"> a. Probe is not equipped with an integral thermocouple b. Low batteries 	<ol style="list-style-type: none"> a. Use a probe that is equipped with a thermocouple b. Replace batteries

WALLACH WA1000B CRYO CONSOLE CLEANING/DISINFECTION INSTRUCTIONS

After each patient use, follow the combined cleaning/disinfecting instructions provided below for the Wallach WA1000B Cryo Console.

1. Do not remove the cryotip from the shaft of the gun
2. Remove disposable shield from cryotip and discard.
3. Instruments should be cleaned as soon as possible. Do not allow blood and debris to dry on the instrument. If cleaning cannot be performed immediately, keep the probe and shaft moist.
4. Prepare an enzymatic, neutral pH cleaner solution according to the Manufacturer's instructions.
5. Wearing appropriate personal protective equipment, hold the cryosurgical gun by the handle with the shaft and probe pointing downwards at an angle. Ensure the probe is still secured to the gun shaft. The probe must remain secured to the gun shaft throughout shaft cleaning.
6. Immerse a clean, lint-free wipe into cleaning solution and thoroughly wring. Thoroughly wipe the black plastic handle including trigger and white insulation probe. Discard the wipe.
7. Immerse a fresh, clean, lint-free wipe into cleaning solution and thoroughly wring. Thoroughly wipe the gauge block, yoke connector and white silicone gas hose. Discard the wipe.
8. Using a soft, nylon brush, scrub difficult-to-access areas such as crevices or textured surfaces.
9. Immerse a fresh, clean, lint-free wipe under warm flowing, utility water and thoroughly wring excess water. Thoroughly wipe the black plastic handle including trigger and white insulation probe for at least 30 seconds. Discard the wipe.
10. Rinse a fresh, clean, lint-free wipe under warm flowing, utility water and thoroughly wring excess water. Thoroughly wipe the gauge block, yoke connector and white silicone gas hose for at least 30 seconds. Discard the wipe.
11. Repeat steps 7 and 8 using a fresh, disposable wipe for each step.
12. Inspect the device for visible soil or debris. If visible soil remains, repeat cleaning steps 4-9 until device is visually clean.
13. Using a fresh moist disinfectant wipe such as Sani-cloth AF3 or CaviWipes, thoroughly wipe the black plastic handle including trigger, shaft (white insulation probe). Discard the wipe.
14. Using a fresh wipe, thoroughly wipe the gauge block, yoke connector and white hose. Discard the wipe.
15. Repeat step 13 with fresh wipes.
16. Discard the wipes and allow the wiped and wetted device to stand for a minimum of three minutes.
17. Use a sterile 70% IPA wipe or a lint-free wipe saturated with sterile 70% IPA, to thoroughly wipe the device.
18. Allow the device to air dry. Once dry, visually inspect the device for any residual soil. Repeat disinfection instruction if visible soil is present.
19. Place the cleaned/disinfected device in an appropriate dry storage area.

WALLACH CRYOSURGICAL TIP CLEANING/STERILIZATION INSTRUCTIONS

After each patient use, follow the combined cleaning/sterilization instructions provided below for the Wallach Cryosurgical Tips.

1. Remove disposable shield from cryotip and discard.
2. Ensure the Cryosurgical Tip is assembled.
3. Prepare an Enzymatic, neutral pH cleaner solution in accordance with Manufacturer's instructions.
4. Grasp the shaft end of assembled device and point tip/probe end at a downward angle.

5. Immerse a clean, lint-free wipe into cleaner solution and thoroughly wring.
6. Beginning just above the shaft joint, thoroughly wipe device with moistened wipe for a minimum of 20 seconds.
7. Use a beaker to pour enzymatic solution, for a minimum of 20 seconds, down the outside of the device, from just above the shaft joint to the tip/probe; do not allow liquid to enter the lumened shaft.
8. After pouring solution, use a clean, soft bristle brush to scrub mated surfaces of the shaft/stem joint and the creases of any welded joints for a minimum of 20 seconds.
9. Rinse the wipe under flowing Utility water and repeat steps 5 through 7.
10. Adjust flowing Utility water to maintain a steady stream at 35-40 °C. Thoroughly rinse device (angled downward), beginning above the shaft joint to the tip/probe, for a minimum of 30 seconds; do not allow liquid to enter the lumened shaft.
11. Pat device dry using a clean, dry lint-free wipe.
12. Holding the device at a downward angle, carefully remove the tip/probe from the shaft.
13. Moisten a lint-free wipe in flowing Utility water and wring thoroughly.
14. Wipe protected area (the stem) of the probe and inside threads of any shaft and stem with moistened wipe. Rinse wipe in flowing Utility water and repeat wiping of these areas for 30 seconds until no visible detergent remains.
15. Pat device dry with a clean low-lint wipe. Visually inspect device for detergent residue.
16. If needed, reconnect tip/probe to shaft and repeat rinse steps 9 through 14 until no visible detergent residue remains.
17. Apply Isopropyl alcohol (70% or greater) to a low-lint wipe (do not saturate). Wipe down the entire device.
18. Dry the device using a clean lint-free wipe.

STERILIZATION

1. Cryotips can be steam sterilized by Gravity or Pre-vacuum using the following procedures:

	Stage / Process	Parameter Set-Points	
Set-up	Cycle Profile	Gravity	Pre-vacuum
Pre-Exposure (Conditioning)	Chamber Purge	1 minute	1 minute
Exposure (Steam Injection)	Sterilization Temperature	132 °C (270 °F)	132 °C (270 °F)
	Sterilization Time	15 minutes	4 minutes
Post-Exposure (Exhaust / Dry)	Dry Time	15 minutes	20 minutes

2. Attach new disposable plastic shield to the cryotip, ensuring that the tabs on the shields are correctly fitted into the grooves on the metal tip.
3. Cryotip is now ready for use.

WARNING: DO NOT USE ETO, E-BEAM OR GAMMA STERILIZATION.

WARRANTY

The WA1500 (Mobile Cart), the WA2000 (Timer/Temperature Indicator Console) and the WA4000 (General Purpose Probe Console) each carry a one-year warranty. **DO NOT ATTEMPT IN-HOUSE REPAIRS AS THEY WILL VOID YOUR WARRANTY.**

SERVICE AND REPAIR

In case of malfunction, immediately turn gas cylinder master valve “OFF”. Also shut off the console switch. Detach console from gas cylinder. Call CooperSurgical Service and Repair at (203) 799-2000 for assistance.

If a repair is needed, carefully sanitize and package the unit in a protective carton. Items that are not sanitized will be returned to the customer freight collect. Include a note stating the problem and ask for a repair estimate. All shipments must be made via prepaid parcel post or U.S. Mail. C.O.D. packages will not be accepted. Address package to the attention of Coopersurgical Service and Repair:

95 Corporate Drive
Trumbull, CT 06611 USA
Phone: 203-799-2000
Fax: 203-799-2002

EXPLANATION OF SYMBOLS

Reorder Number

Serial Number

Consult instructions for use

Caution

Air Flow (Gas Exhaust Port)

System OFF

System ON

R_x Only

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

Authorized Representative in the European Community

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© 2020 CooperSurgical, Inc.