



MEDICAL

**647 + 646**  
**Podiatry Chairs**



# Better position yourself for podiatry procedures.




The integration of new technologies and advancements with in-office procedures are shaping the future of the podiatry practice. Function, efficiency and versatility are key to an effective podiatry room design—and it all starts with the procedure chair.



# Rediscover the most important part of the podiatry chair, **the foot section.**

The **Midmark® 647 Podiatry Chair** with an exclusive foot section design puts control right where it's needed—in the provider's hands. The ergonomic design allows the provider to work from a comfortable position whether seated or standing. Features such as easy-release handles on both sides of the foot section and convenient, integrated controls help make the provider's job a little easier.



## **01 INTEGRATED HAND CONTROLS**

Conveniently located touch controls on both sides of the footrest put operation right at the provider's fingertips.

## **02 ERGONOMIC FOOT SECTION**

Ergonomically positioned release handles allow for one-handed operation. The stainless-steel debris pan is positioned to catch debris while staying out of the way of the provider.

## **03 COVERED GLIDES**

An aseptic design helps simplify cleaning and maintenance by preventing debris and dirt build up. The result is smooth, unobstructed operation.

# Functional. Efficient. Everything within reach.

The Midmark® 647 Podiatry Chair is designed with intelligent features and options to enhance the efficiency of the podiatry procedure. Health systems can custom order chairs to meet the specific needs of their providers and patients while helping improve the overall efficiency of the podiatry procedure workflow.



## ROTATION

Outfitting the chair with optional rotation allows providers to bring the patient and necessary devices to the point of care. No more dragging heavy devices to and from the chair or straining to reach the patient and procedure site.



## ELECTRICAL RECEPTACLES

Optional hospital-grade electrical receptacles and wireless controls free up floor space from power cords, creating a clean look in the procedure room.



## QUICK EXAM®

The standard foot control comes with a Quick Exam button that elevates the chair to the exam position with one touch of the control.



## SEALED FOOT CONTROL

The foot control is sealed against dirt and debris, helping to make cleaning simple and ensure a long product life.



## Patient positioning made easy

### POWERED MULTI-POSITIONING

Our custom-designed motors quickly move your patients into proper position, allowing you to be more productive. Powered movement of the seat, back and tilt gives you the flexibility of positioning patients where you need them for a more ergonomic working position.

### PROGRAMMABILITY

The Midmark 647 can store up to three customizable positions along with a pre-programmed entry position. With the touch of a button, all sections of the chair move simultaneously to quickly and precisely move the patient to the desired position.

### MIDMARK® 646 PODIATRY CHAIR

The 646 is an economical solution for the podiatrist who works exclusively from the seated position. It has the same design features of the 647 with a fixed seat height of 21 inches. The 646 can be customized with any of the available accessories with the exception of the paper roll holder and wireless controls.



Entry Position



Seated Position





Standing Position



Flat Position

# Simple design considerations can improve patient comfort.

Comfortable patients are more relaxed patients. From the patented Premium Comfort System to optional seamless upholstery, the Midmark® 647 is designed to deliver a more relaxed patient, help reduce movement during procedures and improve the care experience.



## 01 LOW HEIGHT

A low 21" seat height and wide, flat transfer surface helps enable patients to transfer onto and off of the procedure chair with minimal assistance which can improve safety in the procedure room.

## 02 CHAIR ARMS

Contoured, upholstered chair arms are designed to give patients a sense of security when accessing the procedure chair. Arms lower with the chair back as well as swivel and lock into multiple positions, even supporting the prone position. The chair arms are removable without the use of tools.

## 03 PREMIUM COMFORT SYSTEM®

The Premium Comfort System uses unique, contoured upholstery to help cradle the patient for a soothing, chair-like feel.





Accessories

# ACCESSORIES

## SWING ARM INSTRUMENT TRAY



10" W x 13 1/2" L

## VISION BLOCK SCREEN



Adjusts vertically  
and horizontally

## WIRELESS FOOT CONTROL



Includes low battery  
audible indicator\*

## PAPER ROLL HOLDER



Holds a 21" x 3 1/2"  
paper roll\*

## PROTECTIVE FOOT COVER



Ships standard with  
upholstery set

## BASE RAIL



Available in 1"

## CASTERS



3 1/2" added to chair height

## RITTER® 253 LED EXAM LIGHT



Available in chair mount,  
wall mount and mobile  
options

\*Not available for the 646



# COLOR OPTIONS

## STANDARD UPHOLSTERY COLORS

	Mist 813
	Stone 814
	Restful Path 858
	Iron Ore 851
	Lunar Gray 845
	Obsidian 857
	Oasis 843
	Soothing Blue 855
	Dream 847
	Mineral 844
	Healing Waters 848
	Shaded Garden 853
	Citrus 849
	Sandy Retreat 850
	Curative Copper 846
	Cranberry 859
	Dark Linen 856
	Latte 860
	Robust Brown 852
	Deep Earth 854

## ULTRAFREE UPHOLSTERY COLORS

	UltraFree Stone 865
	UltraFree Restful Path 868
	UltraFree Obsidian 867
	UltraFree Harbor 861
	UltraFree Wheatgrass 862
	UltraFree Spice 863
	UltraFree Cranberry 869
	UltraFree Dark Linen 866
	UltraFree Latte 870
	UltraFree Branch 864



Length measures 75" and can be increased to 84" by extending foot section.

# SPECIFICATIONS

## 647

Patient Weight Capacity:  
450 lb (204 kg)

Base Height:  
Minimum 21" (53.3 cm)  
Maximum 40" (101.6 cm)

Foot Section Height:  
With foot section retracted  
15.5" to 54" (39.4 cm to 137.2 cm)

With foot section fully extended  
14" to 57" (35.6 cm to 144.8 cm)

Foot Section:  
26" W x 6" H (66 cm x 15.2 cm)

Foot Extension: 9" (22.9 cm)

Debris Tray: 15.1" W x 10.6" H x 1.5" D  
(38.4 cm x 26.9 cm x 3.8 cm)

Back Support Angle Range:  
0° to 80°

Seat Tilt Angle Range: 0° to 30°

Seat Width with Chair Arms Installed:  
37.5" (95.3 cm)

Foot Control:  
14.5" W x 6" H (36.8 cm x 15.2 cm)

Electrical Outlet Rating:<sup>1</sup>  
115 VAC, 3 A

## OPTIONS

**Premium Upholstery Top:**  
Seamless design for easy cleaning (includes a plastic protective foot cover)

**Premium, UltraFree Top:**  
Seamless, PVC-free material gives you the look and feel of leather in an upholstery that is easy to clean and disinfect

**Programmability:**  
Three positions that you customize to your needs

**Electrical Outlets:<sup>1</sup>**  
Duplex, hospital grade (on each side)

**Rotation:<sup>4</sup>**  
Allows 350° of chair rotation and adds 3.5" to the minimum and maximum seat height specifications

## 646

Patient Weight Capacity:  
450 lb (204 kg)

Base Height:  
Fixed height of 21" (53.3 cm)

Foot Section Height:  
With foot section retracted  
18" to 35" (45.7 cm to 88.9 cm)

With foot section fully extended  
15.5" to 37.5" (39.4 cm to 95.3 cm)

Foot Section:  
26" W x 6" H (66 cm x 15.2 cm)

Foot Extension: 9" (22.9 cm)

Debris Tray: 15.1" W x 10.6" H x 1.5" D  
(38.4 cm x 26.9 cm x 3.8 cm)

Back Support Angle Range:  
0° to 80°

Seat Tilt Angle Range: 0° to 30°

Seat Width with Chair Arms Installed:  
37.5" (95.3 cm)

Foot Control:  
14.5" W x 6" H (36.8 cm x 15.2 cm)

Electrical Outlet Rating:  
115 VAC, 3 A

## ACCESSORIES

- Base rail
- Caster base<sup>1</sup>
- Casters<sup>1</sup>
- 3' power cord
- IV pole<sup>3</sup>
- Paper roll holder<sup>4</sup>
- Vision block screen
- Wireless control, linear-style<sup>2,4</sup>
- Replacement plastic protective foot cover
- Swing arm instrument tray, extension top
- Instrument tray, double arm
- Ritter® 253 LED Exam Light<sup>3</sup>

<sup>1</sup> – Not available on rotation models

<sup>2</sup> – Includes foot control (with batteries)

<sup>3</sup> – Requires base rail accessory

<sup>4</sup> – Only available on 647





Designing better care.®



CARB 93129.2 Phase 2 Compliant and TSCA Title VI Compliant

The color examples shown are the best representation of the original material. Actual color may vary slightly. We strongly recommend that you contact Midmark Customer Experience at 1.800.MIDMARK to request a sample before placing your order.

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: [midmark.com/ISO](http://midmark.com/ISO)

For more information, contact your Midmark dealer or call: 1.800.MIDMARK

Outside the USA call: 1.937.526.3662 or visit our website: [midmark.com](http://midmark.com)

© 2023 Midmark Corporation, Versailles, Ohio USA  
Products subject to improvement changes without notice.  
Litho in USA.

007-10362-00 Rev. B1 (7/23)

