

## Platform-Mounted, Heavy-Duty, Parallel Bars



### MODELS 3-2106 & 3-2109

- Parallel bars are mounted to a 1-1/2" thick hardwood plywood base with ultra-durable UV-coated satin finish
- Safety tread at both ends of the platform
- Middle upright with through pin lock for extra support
- Extra-wide 48" width for bariatric patients
- Width adjustable 1-1/2" one-piece, stainless steel handrails
- Easy height adjustment with spring loaded plungers and numbered height strip
- Heavy duty, square telescoping steel uprights
- Second height knob removes side to side play
- End bumpers on all handrails
- 500 lbs. load capacity under normal physical therapy use

**IMPORTANT NOTE:** *Intended for physical therapy use only, not for gymnastics use.*

### 3-2106 SPECIFICATIONS


10' L x 20"-34" W range x 26"-39" H range

### 3-2109 SPECIFICATIONS

7' L x 20"-34" W range x 26"-39" H range



California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1