

## Panel Leg Recovery Couch with Full Shelf

### MODEL 3710

- Easy-clean all laminate construction
- Full shelf for storage
- 2" firm foam padding
- 250 lbs. load capacity under normal use



**3710-16 Adjustable pillow wedge headrest**

72"L x 27" W x 20"H



**3710-15 w/Flat Foam Adjustable Headrest**

72"L x 27" W x 20"H



California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2




**3710-10w/Non-Adjustable Pillow Wedge**

72"L x 27" W x 20"H

### Laminate Colors



Maple



Dark Cherry



Gray



Slate Gray

### Standard Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black