

KD Panel Leg Series Couch with Drawers

MODEL 3913-27

- All laminate construction
- Sleek, easy-clean, panel leg design
- Top is easily replaceable
- Two drawer storage
- Adjustable backrest with 6 settings of approximately: 27°, 36°, 43°, 50°, 55° & 60°
- Scratch resistant topcoat finish
- Ships flat for easy handling and extra savings
- Fast and easy assembly
- Seamless rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- 300 lbs. load capacity under normal use (181.43 kg)



SPECIFICATIONS

3913-27

72" L x 27" W x 21.5" H



California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

Options



034– Paper Dispenser & Holder Combo



040– Paper Cutter



030– Paper Dispenser



020– Small Pillow



060– Large Pillow



055– Drawer Locks


Laminate Colors



Gray



Maple



Dark Cherry

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black