


## Wedge



### MODEL 56

- Firm high density urethane foam core
- 18° incline
- Knit-backed, sewn, vinyl cover with hook and loop closure
- Easy-clean surface

### MODEL 56 SPECIFICATIONS

24" W x 18" D x 6" H at tall end



**NOTE:** If no vinyl color is selected, Royal Blue will be used.

California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

### Standard Upholstery



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black