

## Imagination Series, Rainforest Follies, Pediatric Scale Table

### MODEL 7832

#### TABLE FEATURES

- Clinton original illustrations
- Easy to clean laminate surfaces
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Graphics are on one long side and two short ends. Back is slate gray. Only manufactured as pictured.
- Comfortable work height of 36"
- Top is rimmed on three sides for safety
- Top's front and back rim edges are beveled and the corners rounded for safety
- Easy-clean, replaceable center pad
- Built-in infantometer with folding ends for measuring height
- Easy access to weight, measurement and storage
- Choose to access concealed paper dispenser from either end
- Recessed front toe kick
- Easy-clean all laminate base
- Ample 2 door and 1 drawer storage
- Soft-close, door hinges



#### MODEL 7832 SPECIFICATIONS

43"L x 24"W x 38.5"H

#### SCALE FEATURES

- Dependable, accurate Doran® scale
- Easy-read digital display
- Auto calibration on start-up
- Auto shut-off
- Electric or battery operation (6 AA batteries required)
- Weights in increments of lbs., ozs. or kilos
- Minimum registered weight is 5 lbs.—Maximum weight for scale is 40 lbs.

California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1



### Options



78-P Parrot Scale Pal



066 Door & Inside Latch Combo




055 Drawer Locks

### Upholstery



Gunmetal

## Features



### Infantometer

Built-in infantometer with folding ends for measuring height



### Paper Dispenser

Choose to access concealed paper dispenser from either end



### Rimmed Top with Beveled Edges

Top is rimmed on three sides for safety and the front and back rim edges are beveled and the corners rounded for safety



### Padded Top

Easy-clean, padded, upholstered top for comfort



### Scale Head

Easy-read digital display