


Imagination Series, Alley Cats & Dogs, Treatment Table



MODEL 7937 (Adjustable Backrest)

MODEL 7937-1 (Flat Top)

- No upcharge for adjustable backrest
- Clinton original themes
- Easy to clean laminate surfaces
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- 2" of firm foam padding
- 4 door storage with 2 adjustable shelves
- Scratch and mar resistant surfaces
- Chip resistant PVC edges
- Premium, Gunmetal, stain-resistant, woven, knit-backed vinyl upholstery
- Soft-close, door hinges
- Graphics are on one long side with slate gray ends and back
- Top features seamless rounded outside corners
- Paper dispenser included
- Available with or without backrest

SPECIFICATIONS (BOTH MODELS)

68"L x 27"W x 33"H

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.


1



Options




066 Door & Inside Latch Combo




055 Door Locks



040 Paper Cutter



20 Small Pillow



60 Large Pillow

Upholstery



Gunmetal