

## Style Line, Cabinet Style, Treatment Table with 4 Doors



### MODEL 9074

- No upcharge for adjustable backrest
- Features 4 doors, adjustable backrest and 2 storage compartments each with an adjustable shelf
- Easy-clean laminate surfaces
- Adjustable doors with concealed soft-close hinges
- Ships assembled
- Seamless rounded corner top
- Heavyweight, knit-backed upholstery
- 2" firm foam padding
- 425 lbs. load capacity under normal use

**NOTE:** When ordering 066 locks, at least 2 locks must be ordered to lock each section (front and back). Locks are not intended to secure narcotics or controlled items.



### SPECIFICATIONS

9074-27	72" L x 27" W x 31" H
9074-30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

### Options



### Laminate Colors



Dark Cherry



Maple



Gray

### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black