

ZOLL AED 3

The AED That Leads the Way



Enhanced Real CPR Help® — our proven real-time CPR feedback technology now features a full-color display with vivid rescue images, a CPR cycle timer, and a large color bar gauge that shows CPR compression depth.

Innovative Design and Enhanced Features

The ZOLL AED 3® defibrillator builds on our reputation for providing rescuers with unrivaled rescue support and a low total cost of ownership. The ZOLL AED 3 embodies all of the same great qualities of the AED Plus® defibrillator, while offering new, advanced features, including:

- **Real CPR Help** — Provides integrated, real-time feedback on CPR compression rate and depth.
- **Integrated pediatric rescue** — Universal CPR Uni-padz® electrodes can be used to treat both adult and pediatric victims of sudden cardiac arrest (SCA) by simply activating child mode.
- **WiFi connectivity** — Automatically reports AED readiness.
- **Long-life consumables** — Five-year battery and pads reduce AED device maintenance requirements and the total cost of ownership of your ZOLL AED 3.

ZOLL AED 3 Package:

Package includes ZOLL AED 3, CPR Uni-padz (for adult and pediatric use), ZOLL AED 3 battery pack, and manuals. Free on-site in-service with every purchase.

	Part #	Sale Price
Semi-automatic ZOLL AED 3 Package		
	Part #	Sale Price
Automatic ZOLL AED 3 Package		



350,000

Out-of-hospital cardiac
arrests in 2016¹



10%


Decrease in survival for every
minute of delayed defibrillation²



5%

Chance of survival if a victim
collapses and nothing is done
until emergency services arrive³

Product Listings and Accessories:



CPR Uni-padz

(8900-000260):

Truly universal pads. Use the same pads in an adult or child rescue by simply repositioning the pads and selecting adult or child rescue mode. Five-year installed life. Provides full CPR feedback for adult patients.



ZOLL AED 3 Battery Pack

(8000-000696):

In stand-by mode, the ZOLL AED 3 disposable lithium manganese dioxide battery has an installed life of up to five years.



Standard Carry Case

(8000-001252):

In addition to the ZOLL AED 3 with preconnected CPR Uni-padz, the Standard Carry Case can also hold a backup set of CPR Uni-padz.



Surface-mounted Cabinet

(8000-001256):

Holds a ZOLL AED 3 plus a backup set of electrodes and a backup battery.

*Additional cabinet versions and wall mounts available



Large Hard-shell Carry Case

(8000-001254):

Protect and carry everything needed for a vigorous rescue with the large hard-shell carry case, including a backup set of CPR Uni-padz and a backup battery.



AED 3 Trainer

(8028-000001-01):

The trainer will allow students to learn proper rate and depth of CPR chest compression in a training setting. Instructors can set up various rescue scenarios.



Want to know more?

Please contact your ZOLL Medical Distributor:

¹American Heart Association website. http://cpr.heart.org/AHA/ECC/CPRandECC/General/UCM_477263_Cardiac-Arrest-Statistics.jsp. Accessed 21 July 2020.

²Weisfeldt ML, et al. *J Am Coll Cardiol*. 2010;55(16):1713–20.

³*Circulation*. 2010;122:S720–S728, <https://doi.org/10.1161/CIRCULATIONAHA.110.970970>.