

## Gray Frame, Exam Room Chair/No Arms

### MODEL C-40G

- Gray powder coated frame finish
- Heavy duty 1-1/2" diameter legs
- Solid plywood seat and back
- Wall saver guard
- Contoured seat and back
- 2" thick, firm foam padded seat and back
- Adjustable leg levelers
- Premium, knit-backed upholstery
- 300 lbs. load capacity under normal use



### SPECIFICATIONS

Seat	18" W x 18" D
Seat Height	18"

California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals including 1, 3-Butadiene which is known to the State of California to cause cancer and reproductive harm. For more information go to: [www.P65Warning.ca.gov](http://www.P65Warning.ca.gov).

7

### Upholstery Colors

