

Flat Top, Classic Series, Straight Line Treatment Table



MODEL 100

- No upcharge for flat top
- Ultra-durable, PVC-laminated, hardwood legs
- H-brace construction
- Legs have pre-inserted carriage bolts for fast easy assembly
- Seamless, rounded corner top
- 2" of firm foam padding (5 cm)
- Available in Natural and Dark Cherry finish
- Choice of upholstery colors
- Fast, easy assembly with basic tools
- 400 lbs. load capacity under normal use (181 kg)

SPECIFICATIONS

100-24	72" L x 24" W x 31" H
100-27	72" L x 27" W x 31" H
100-30	72" L x 30" W x 31" H
100-38	78" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Options



20- small Pillow



034- Paper Dispenser & Cutter Combo



040- Paper Cutter




030- Paper Dispenser



005-Special Height



006-Special Length



009-Nose Hole

Finishes



Natural



Dark Cherry

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black