

Classic Series, Straight Line Treatment Table with Drawers




MODEL 1013

- No upcharge for adjustable backrest
- Ultra-durable, PVC-laminated, hardwood legs
- Laminate shelf, two easy clean, plastic drawers
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Drawers can be positioned at either end, and from either side
- Scratch resistant topcoat finish
- Fast and easy assembly
- Seamless rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- 400 lbs. load capacity under normal use (181 kg)

SPECIFICATIONS

1013-27	72" L x 27" W x 31" H
1013-30	72" L x 30" W x 31" H
1013-38	78" L x 30" W x 31" H


Prop 65 Warning



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Options



20- small Pillow



034- Paper Dispenser & Cutter Combo



040- Paper Cutter




030- Paper Dispenser



005-Special Height



006-Special Length



009-Nose Hole




055-Drawer Locks

Finishes



Natural



Dark Cherry

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black