

Classic Series, Straight Line Treatment Table with Shelving



MODEL 1030

- No upcharge for adjustable backrest
- Ultra-durable, PVC-laminated, hardwood legs
- Laminate shelf and storage compartments
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Scratch resistant topcoat finish
- Fast and easy assembly
- Seamless rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- 400 lbs. load capacity under normal use (181 kg)

SPECIFICATIONS


1030-27	72" L x 27" W x 31" H
1030-30	72" L x 30" W x 31" H
1030-38	78" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Options




20–small Pillow



034– Paper Dispenser & Cutter Combo



040– Paper Cutter



030– Paper Dispenser



005–Special Height



006–Special Length



009–Nose Hole

Finishes



Natural



Dark Cherry

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black