

Chrome Base Stool

MODEL 2102

- 2" rubber wheel, ball bearing casters
- solid steel machine screw adjustment
- All models have 19.5" diameter bases
- Cover cap conceals top welds
- Seat is mounted to a 6" x 6" steel base
- Premium, stain resistant, woven, knit-backed vinyl upholstery
- 250 lbs. load capacity under normal use

SPECIFICATIONS

Seat Diameter	14"
Base Diameter	19.5"
Height Range	18.5"-25.5"



California residents please be advised of the following, pursuant to Proposition 65:

! **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Options



24 Optional Rubber Tips



23 Optional Glides C

Upholstery Colors

