

## 5-Leg Pneumatic Stool

### MODEL 2156

- Premium 23" cast aluminum base with satin top and black accents
- 5-leg design for greater stability
- Full size diameter top and base for comfort and stability
- 4" thick seat for comfort
- Chrome hooded, dual wheel, nylon casters
- Large, easy to locate, D-shaped height lever
- Premium, stain resistant, woven, knit-backed upholstery
- Knock-down design saves on shipping
- Easy, no-tools assembly
- Also available with backrest
- 250 lbs. load capacity under normal use



### SPECIFICATIONS

Seat Diameter	16"
Base Diameter	23"
Height Range	18.5"–23.5"

California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

### Options



23 Optional Glides B

25 Soft-Roll Casters

27 Scuff Protectors

22 Optional Foot Ring

26 Extra Tall Pneumatic Cylinder

### Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black