

## Alpha Series Treatment Table with H-Brace



### MODEL 3010

- No upcharge fo adjustable backrest
- Gray, powder-coated, all-steel frame, and adjustable backrest
- Easy bolt-through assembly
- Paper dispenser included
- Adjustable leg levelers
- Seamless, rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- 450 lbs. load capacity under normal use (204 kg)

### SPECIFICATIONS

3010-27	72" L x 27" W x 31" H
3010-30	72" L x 30" W x 31" H

**California residents please be advised of the following, pursuant to Proposition 65:**


**!** **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

**2**

### Options



020 Small Pillow



060 Large Pillow



009 Nosehole



040 Paper Cutter

### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black