

## Flat Top, Alpha Series, Treatment Table with Drawers



### MODEL 3130

- No upcharge for flat top
- Gray, powder-coated, all-steel frame, gray laminate shelf, and two, easy-clean, plastic drawers with gray, drawer fronts
- Adjustable height levelers
- Pre-threaded hardware for easy assembly
- Paper dispenser included
- 2" of firm foam padding (5 cm)
- Seamless, rounded corner, flat top
- Fast, easy assembly with basic tools
- 450 lbs. load capacity under normal use (204 kg)

### SPECIFICATIONS

3130-27	72" L x 27" W x 31" H
3130-30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:


**!** **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1


### Options



020 Small Pillow



060 Large Pillow



009 Nosehole



040 Paper Cutter



055 Drawer Locks

### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black