

600 Series, Metal, Treatment Table



MODEL 6198

- Features gray, powder-coated, all-steel frame
- Plywood top and welded steel frame
- 6-leg design with bolted-through assembly
- Leg levelers
- Gray, laminate, full length shelf with H-braces
- Seamless, rounded corner, flat top
- 2" firm foam padding (5 cm)
- 600 lbs. capacity under normal use (272 kg)

SPECIFICATIONS

6198	72" L x 30" W x 31" H
------	-----------------------

California residents please be advised of the following, pursuant to Proposition 65:

! **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Options



020 Small Pillow



060 Large Pillow



009 Nosehole



040 Paper Cutter

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black