

H Series, Padded, Blood Drawing Chair with Padded Arms

MODEL 6310

- Padded seat and backrest
- Seat is slanted back 6° for comfort and stability
- Height and depth adjustable, padded, flip arm
- Flip arm positions in place when in use, then flips out of the way for easy exit
- Height-adjustable armrest
- Padded, elbow rest
- 1-1/4" square, heavy duty, all-welded, tubular steel, neutral gray powder-coated seat frame (3.18 cm)
- Polished chrome base with black bumper
- Large foot pedal
- Easy-clean, durable vinyl in a selection of colors
- 300 lbs. load capacity under normal use (136 kg)

IMPORTANT NOTE: Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

ORDERING NOTE: MUST add arm suffix number to model number when ordering with 665 arm option, otherwise the chair will be shipped as shown.

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.



MODEL 6310 SPECIFICATIONS

Seat	21" W x 19.5" D
Overall	30.75" W x 30" D
Height Range	20.25" — 27.5"
Arm Adjustment	7"



California residents please be advised of the following, pursuant to Proposition 65:

1
WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

OPTIONS



665 ClintonClean™, Armrests and Flip Arm(s)

Standard Upholstery Colors

