

H Series, E-Z-Clean, Blood Drawing Chair with ClintonClean™ Flip Arm and Drawer

MODEL 6340-P

- E-Z-Clean, one-piece, plastic seat
- Seat is slanted back 6° for comfort and stability
- Height and depth adjustable, ClintonClean™ flip arm
- Flip arm positions in place when in use, then flips out of the way for easy exit
- ClintonClean™ elbow rest
- ClintonClean™ drawer top (*NOTE: Due to the 6° slant, objects should not be placed on the drawer top*)
- Height-adjustable, one-piece, all-plastic drawer
- 1-1/4" square, heavy duty, all-welded, tubular steel, neutral gray powder-coated frame (3.18 cm)
- Polished chrome base with black bumper
- Large foot pedal
- Height adjustable armrest
- 300 lbs. load capacity under normal use (136 kg)


NOTE: If padded upholstered arms are selected, use color chart to choose color.

IMPORTANT NOTE: Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

ORDERING NOTE: MUST add arm suffix number to model number when ordering with 600 arm option, otherwise the chair will be shipped as shown.

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.

California residents please be advised of the following, pursuant to Proposition 65:

1
 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.



MODEL 6340-P SPECIFICATIONS

Seat	18" W x 17" D
Overall	38" W x 30" D
Height Range	20.25" — 27.5"
Arm Adjustment	7"

600 OPTION



Padded, Armrest and Flip Arm



Standard Upholstery Colors for Optional Arms

