

Power Series, Extra-Wide, Blood Drawing Chair with Padded Flip Arm & Drawer

MODEL 6352

- Padded seat and backrest
- Height and depth adjustable, padded flip arm
- Flip arm positions in place when in use, then flips out of the way for easy exit
- Padded, elbow rest
- ClintonClean™ drawer top
- Height-adjustable, one-piece, all-plastic drawer
- 1-1/4" square, heavy duty, all-welded, tubular steel, neutral gray powder-coated frame (3.18 cm)
- Contemporary brushed aluminum motor casing
- Hands-free, foot pedal controller
- Easy-clean, durable vinyl in a selection of colors
- 375 lbs. load capacity under normal use (170 kg)

IMPORTANT NOTE: Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

ORDERING NOTE: MUST add arm suffix number to model number when ordering with 665 arm option, otherwise the chair will be shipped as shown.

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.



MODEL 6352 SPECIFICATIONS

Seat	30" W x 18" D
Overall	48" W x 30" D
Height Range	21.5" — 29.5"
Arm Adjustment	8"



California residents please be advised of the following, pursuant to Proposition 65:

1
WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

OPTIONS



665 ClintonClean™, Armrest and Flip Arm(s)

Standard Upholstery Colors

