

Power Series, Bariatric, Blood Drawing Chair with Padded Arms

MODEL 6361

- Padded seat and backrest
- Dual, height and depth adjustable, padded, flip arms
- Flip arm(s) positions in place when in use, then flips out of the way for easy exit
- Dual, padded, elbow rests
- 1-1/4" square, heavy duty, all-welded, tubular steel, neutral gray powder-coated frame (3.18 cm)
- Contemporary brushed aluminum motor casings
- Hands-free, foot pedal controller
- Easy-clean, durable vinyl in a selection of colors
- 850 lbs. combined load capacity under normal use (385.55 kg)

IMPORTANT NOTE: Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

ORDERING NOTE: MUST add arm suffix number to model number when ordering with 665 arm option, otherwise the chair will be shipped as shown.

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.



MODEL 6361 SPECIFICATIONS

Seat	35" W x 20" D
Overall	45" W x 32" D
Height Range	21.5" — 29.5"
Arm Adjustment	8"



California residents please be advised of the following, pursuant to Proposition 65:

1
WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

OPTIONS



665 ClintonClean™, Armrests and Flip Arm(s)

Standard Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black