

## Lab X Series, Extra-Wide, Blood Drawing Chair with Padded Arms

### MODEL 66000

- Padded seat and backrest
- Height and depth adjustable, dual, padded, flip arms
- Flip arm(s) positions in place when in use, then flips out of the way for easy exit
- Dual ClintonClean elbow rests
- ClintonClean™ elbow rests
- A full 35" wide seat to accommodate a larger patient, or parent and child (89 cm)
- 1-1/4" square, heavy duty, all-welded, tubular steel frame (3.18 cm)
- Comfortable 2" padding throughout (5 cm)
- Neutral Gray, powder-coated frame
- Adjustable foot levelers
- Easy-clean, durable vinyl in a selection of colors
- ClintonClean™ arms option available
- 400 lbs. combined load capacity under normal use (181.4 kg)




**IMPORTANT NOTE: Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.**

**ORDERING NOTE: MUST add arm suffix number to model number when ordering any model with arm option, otherwise chair will be shipped as shown.**

**All blood drawing products are non-returnable unless they are unused in unopened, original packaging.**



California residents please be advised of the following, pursuant to Proposition 65:

**1**  
 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

### MODEL 66000 SPECIFICATIONS

|                  |                       |
|------------------|-----------------------|
| Seat             | 35" W x 18" D x 20" H |
| Overall          | 46" W x 26" D x 36" H |
| Arm Height Range | 27" — 33"             |
| Arm Adjustment   | 6"                    |

### OPTIONS



665 ClintonClean™, Armrests and Flip Arm(s)

### Standard Upholstery Colors

