

Lab X Series, Extra-Wide, Blood Chair with Padded Flip Arm and Drawer

MODEL 66002

- Padded seat and backrest
- Height and depth adjustable, padded flip arm
- Flip arm(s) positions in place when in use, then flips out of the way for easy exit
- ClintonClean™ elbow rest
- ClintonClean™ drawer top
- Height-adjustable, one-piece, all-plastic drawer
- A full 35" wide seat to accommodate a larger patient, or parent and child (89 cm)
- 1-1/4" square, heavy duty, all-welded, tubular steel frame (3.18 cm)
- Comfortable 2" padding throughout (5 cm)
- Neutral Gray, powder-coated frame
- Easy-clean, durable vinyl in a selection of colors
- ClintonClean™ arms option available
- 400 lbs. combined load capacity under normal use (181.4 kg)



IMPORTANT NOTE: Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

ORDERING NOTE: MUST add arm suffix number to model number when ordering any model with arm option, otherwise chair will be shipped as shown.

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.



California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

MODEL 66002 SPECIFICATIONS

Seat	35" W x 18" D x 20" H
Overall	53" W x 26" D x 36" H
Arm Height Range	27" — 33"
Arm Adjustment	6"

OPTIONS



665 ClintonClean™, Armrest and Flip Arm(s)

Standard Upholstery Colors

