

Lab X Series, Blood Drawing Chair with Padded Flip Arm and Drawer

MODEL 66020

- Contoured padded seat and padded, arched backrest
- Height and depth adjustable, padded, flip arm
- Flip arm positions in place when in use, then flips out of the way for easy exit
- ClintonClean™ elbow rest
- ClintonClean™ drawer top
- Height-adjustable, one-piece, all-plastic drawer
- 1-1/4" square, heavy duty, all-welded, tubular steel, frame (3.18 cm)
- Comfortable, 2" padding throughout (5 cm)
- Neutral Gray, powder-coated frame
- Adjustable foot levelers
- Easy-clean, durable vinyl in a selection of colors
- ClintonClean™ arms option available
- 400 lbs. load capacity under normal use (181.4 kg)



IMPORTANT NOTE: Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

ORDERING NOTE: MUST add arm suffix number to model number when ordering with 665 arm option, otherwise the chair will be shipped as shown.

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.

MODEL 66020 SPECIFICATIONS

Seat	17.5" W x 18" D x 20" H
Overall	40" W x 27" D x 36"H
Arm Height Range	27" — 33"
Arm Adjustment	6"



California residents please be advised of the following, pursuant to Proposition 65:

1
WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

OPTIONS



665 ClintonClean™, Armrest and Flip Arm(s)

Standard Upholstery Colors

