

## Lab X Series, Extra-Tall, Blood Drawing Chair w/ Padded Flip Arm and Drawer

### MODEL 66022

- Padded seat and padded, arched backrest
- Height and depth adjustable, padded, flip arm
- Flip arm positions in place when in use, then flips out of the way for easy exit
- ClintonClean™ elbow rest
- ClintonClean™ drawer top
- Height-adjustable, one-piece, all-plastic drawer
- Extra height reduces bending for attendant
- 1-1/4" square, heavy duty, all-welded, tubular steel, frame (3.18 cm)
- Comfortable, 2" padding throughout (5 cm)
- Neutral Gray, powder-coated frame
- Easy-clean, durable vinyl in a selection of colors
- ClintonClean™ arms option available
- 400 lbs. load capacity under normal use (181.4 kg)

**IMPORTANT NOTE:** Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

**ORDERING NOTE:** MUST add arm suffix number to model number when ordering with 665 arm option, otherwise the chair will be shipped as shown.

**All blood drawing products are non-returnable unless they are unused in unopened, original packaging.**



MODEL 66022 SPECIFICATIONS

Seat	17.5" W x 18" D x 26" H
Overall	40" W x 27" D x 42"H
Arm Height Range	33" — 39"
Arm Adjustment	6"



California residents please be advised of the following, pursuant to Proposition 65:

**1**  
**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

### OPTIONS



665 ClintonClean™, Armrest and Flip Arm(s)

### Standard Upholstery Colors

