

Lab X Series, Heavy-Duty, Extra-Tall, Blood Drawing Chair with Padded Arms

MODEL 66099B

- Padded, larger seat features extra width and depth for added comfort
- Comfortable padded arched backrest
- Height and depth adjustable, dual, padded, flip arms
- Flip arm(s) positions in place when in use, then flips out of the way for easy exit
- Dual, ClintonClean™ elbow rests
- A full 35" wide seat (89 cm)
- 1-1/4" square, heavy duty, all-welded, tubular steel frame (3.18 cm)
- Comfortable 2" padding throughout (5 cm)
- Neutral Gray powder-coated frame
- Easy-clean, durable vinyl in a selection of colors
- ClintonClean™ arms option available
- 700 lbs. combined load capacity under normal use (317.8 kg)



IMPORTANT NOTE: Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

ORDERING NOTE: MUST add arm suffix number to model number when ordering any model with arm option, otherwise chair will be shipped as shown.

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.

MODEL 66099B SPECIFICATIONS

Seat	35" W x 21" D x 26" H
Overall	46" W x 29" D x 42.5"H
Arm Height Range	33" — 39"
Arm Adjustment	6"

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

OPTIONS



665 ClintonClean™, Armrests and Flip Arm(s)



Standard Upholstery Colors

