

Power 600, Shrouded, Extra Wide, Heavy-Duty, Straight Top Power Table

MODEL 84108

- Extra-wide model with choice of two widths
- Features easy-clean, one-piece top
- Easy-clean ABS plastic enclosed lift
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Rounded corner, upholstered top
- 2" (5 cm) ultra-firm padding
- Customize with options
- Made in the USA by Clinton Industries
- 600 lbs. (272 kg) load capacity under normal use

POWER 600




84108-34 SPECIFICATIONS

Overall 72" L x 34" W x 21"-34" H

84108-40 Extra-Wide 40" width model

California residents please be advised of the following, pursuant to Proposition 65:

1
 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.



Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

OPTIONS

Safety Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 safety straps per set (92" x 2")

080



Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

082



Hand-Held Controller

- Alternative to foot control
- Stores easily under top when not in use

084



Locking Casters

- 3" Heavy duty, rubber wheel, locking, swivel casters on all four corners makes it easy to reposition the table and then lock it securely in place.

087



Paper Dispenser

- Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

030



Paper Cutter

- Clear plastic with elastic straps and buckle

040



Paper Dispenser/Cutter Combo

- Paper dispenser and paper cutter ordered as one option

034