

Power 600, Shrouded, Extra Wide, Heavy-Duty, Power Table with Adjustable Backrest and Drop Section

MODEL 84438

- Extra-wide model with choice of two widths
- Easy-clean ABS plastic enclosed lift
- Pneumatic backrest is adjustable from either side and raises 80 degrees from flat
- Drop Section is adjustable from either side and lowers up to 80 degrees from flat
- Standard drop section has nose hole, but it may be ordered without a nose hole (Please specify when ordering)
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Rounded corner, upholstered top
- 2" ultra-firm padding
- Customize with options
- Made in the USA by Clinton Industries
- 600 lbs. (272 kg) load capacity under normal use

POWER 600



84438-34 SPECIFICATIONS

Overall 76" L x 34" W x 21"-34" H

Backrest 38" L x 34" W

Drop Section 13" L x 34" W

84438-40 Extra-Wide 40" width model

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

OPTIONS

Safety Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 safety straps per set (92" x 2")

080



Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

082



Hand-Held Controller

- Alternative to foot control
- Stores easily under top when not in use

084



Locking Casters

- 3" Heavy duty, rubber wheel, locking, swivel casters on all four corners makes it easy to reposition the table and then lock it securely in place.

087



Paper Dispenser

- Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

030



Paper Cutter

- Clear plastic with elastic straps and buckle

040



Paper Dispenser/Cutter Combo

- Paper dispenser and paper cutter ordered as one option

034