

Summit Doppler™

LifeDop®

A Complete Family of OB Dopplers

The Summit Doppler™ LifeDop® family of Dopplers offers superior sensitivity, high-fidelity sound, a recording option, and state-of-the-art technology.



L150 Series – The LifeDop150 is a non-display handheld Doppler with superior sound quality that accepts all non-directional LifeDop probes. Rechargeable batteries and audio recording options available on specific models.

L250 Series – The LifeDop 250 is a full-featured handheld device with a large display that indicates heart rate, battery status, and signal strength. Rechargeable batteries and audio recording options are available on specific models. The LifeDop 250 accepts all LifeDop non-directional probes.

L350 Series – The LifeDop 350 is a convenient tabletop Doppler offering excellent sound quality with its high-fidelity speaker. It features a large heart rate display, rechargeable batteries, and accepts all LifeDop non-directional probes. Its enhanced features make it more efficient for detecting fetal heartbeats.

LifeDop Probes – Sensitive, Durable, Interchangeable

All LifeDop obstetrical probes have built-in electronics that provide excellent sensitivity, complete interchangeability, and maximum durability. Choose the best probe for your application with the freedom to add additional probes in the future, all on the same unit.

2 MHz – This wide-angle beam is optimized for deep penetration for late term pregnancy and larger patients.

2 MHz Waterproof – Designed for water-assisted labor and delivery applications.

3 MHz – A general-purpose obstetrical probe for early fetal heartbeat detection through delivery. Wide-angle beam decreases time spent “hunting” for the heartbeat. Once fetal heart rate is detected, the sound is clearer than ever.

Summit Doppler™

Probes – All of our Dopplers have interchangeable probe capabilities so you can choose the best probe for your needs



LifeDop® 150 Series



LifeDop® 250 Series

Main Units	Interchangeable Probes	Superior Sound Quality	Rechargeable Batteries	Graphic LCD Display	Built-in Audio Recorder
L150 Non-Display Doppler	X	X			
L150R Non-Display Doppler with Recharger	X	X	X		
L150A Non-Display Doppler with Audio Recorder	X	X			X
C250CB Display Doppler, 2 & 3 MHz Probes	X	X		X	
L250 Display Doppler	X	X		X	
L250R Display Doppler with Recharger	X	X	X	X	
L250AR Display Doppler with Recharger and Audio Recorder	X	X	X	X	X
L350R Table-Top Display with Recharger and High-Fidelity Speaker	X	X	X	X	

Item Number	Description
SD2	2 MHz Obstetrical Probe – Optimized for deep penetration
SDW	2 MHz Obstetrical Waterproof Probe – Optimized for deep penetration and water applications
SD3	3 MHz Obstetrical Probe – For early fetal heartbeat detection and overall use



LifeDop® 350 Series

Stand is an optional accessory

CooperSurgical

800.243.2974 | 203.601.5200
www.coopersurgical.com/distribution
 95 Corporate Drive, Trumbull, CT 06611