

DETECTO®

6868 | BARIATRIC FLIP SEAT SCALE

Digital Multi-Functional Scale with Integrated RS232 Interface



- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- Unique flip-up seat for standing or seated weighing
- Extra-large 32 in x 34 in (81 cm x 86 cm) platform
- Padded stainless steel wrap-around handrails
- Wheels and balanced design for mobility



*OIML EC Class III model
also available (6868KGEU)*

6868

BARIATRIC FLIP SEAT SCALE

	6868	6868KGEU
Capacity	1,000 lb x 0.2 lb / 450 kg x 0.1 kg	300 kg x 0.1 kg or 350 kg x 0.2 kg (OIML EC Class III)
DIMENSIONS	Overall Dimensions	32.3 in W x 39.3 in D x 44.0 in H / 81.9 cm W x 99.7 cm D x 111.8 cm H
	Platform Dimensions	32.0 in W x 33.7 in D x 2.4 in H / 81.3 cm W x 85.7 cm D x 6.0 cm H
Construction	Aluminum and stainless steel	
Power	6 AA batteries (not included) or optional AC adapter	6 C batteries (not included) or optional AC adapter
	Display	Extra-large, 0.75 in/ 19 mm-high LCD digits
INDICATOR	Keypad	Membrane type with 5 keys and 4 directional arrows
	Buttons	Zero, Units, Lock/Release, On/Off, BMI/Enter, Print, Arrow keys
	Connectivity	RS232 serial port
WEIGHT	Net Weight	80 lb / 36 kg
	Shipping Wt.	178 lb / 81 kg



The 6868 clinical scale features an easy-to-use keypad and extra-large display for measurement readouts.



The scale's RS232 serial port provides EMR/EHR connectivity for efficient patient record-keeping.



The 6868 bariatric scale uses 6 AA batteries (not included) or an optional AC adapter (model 6800-1045).

DETECTO's model 6868 clinical scale offers a variety of usages due to its unique product design and high capacity: standing, seated, and bariatric weighing. The wraparound stainless steel handrails with padded grips provide comfortable support for obese and unsteady patients. The 6868 is completely portable and can be transported easily due to its battery power and integral wheels (6 AA batteries not included).

- Body Mass Index calculation
- Lock/Release key for viewing weight after the patient steps off the scale
- Automatic Shutoff and Sleep modes conserve battery life
- Sturdy, long-lasting quality construction
- Many other bariatric scales also available from DETECTO
- Fits through doorways up to 32 inch (81 cm) wide

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.

SOLD BY:

DETECTO®

102 East Daugherty St., Webb City, MO 64870 USA
Toll-Free: (800) 641-2008 E-mail: detecto@cardet.com
www.Detecto.com