

FDA Cleared

enbio.

Unwrapped
FAST / 273°F


7 min

Wrapped
Vacuum / 273°F

15 min

Wrapped
Vacuum / 250°F

45 min



Enbio S Autoclave Technical Data




Vacuum
Pouches+
drying

Capacity



x8



Load space
11×7×1.3 in

Max. packages amount (circa) **8**
Max. recommended batch weight **1.1 lb**
Process chamber volume **2.7 l**



Power supply

Enbio S works with any typical electrical outlet.



Data archiving

All process data are automatically written to a USB disk.



Water supply

Enbio uses distilled water from external tank.



Process reports

Take full control of your sterilization process with my.enbio app.

Dimensions



22 in



9.8 in

6.3 in

Pressure & temperature

Working pressure **2.1 Bar/30.5 psi**
Max. pressure **2.45 Bar/33.4 psi**
Max. process temperature **137°C/278.6°F**

Electrical

Power supply **110-120V/60Hz**
Nominal power **3.25 kW**
Max. power consumption **15 A**

Misc

Ingress protection **IP20**
Data archiving **USB**
Water conductivity* **< 15 µS/cm**

* - demineralized water only

Enbio Group AG
Eichengasse 3, 4702 Oensingen
Switzerland

US Distribution Center: Enbio Corp.
1075 Highway 124, Building 300, Suite 300
Hoschton, GA 30548

enbio.com