

**Health o meter®**  
P r o f e s s i o n a l

10 Year Limited  
Warranty



### Advantages

- ✓ **Powder Coated Fabrication**  
Powder coating inside and out protects it from rust, extending the life of the scale
- ✓ **Exclusive 10 Year Limited Warranty**  
Health o meter® Professional's industry leading warranty protects you against defects of materials or workmanship

**Health o meter®**  
 P r o f e s s i o n a l



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
400KL	58.875" x 20.25" x 6.375" / 4.40	28.75	892076002121

**Specifications**

- Capacity: 390 lb / 180 kg
- Graduation: 1/4 lb / 100 g
- Platform Size: 10-1/2" (w) x 14" (d) x 3-1/4" (h), 267 mm (w) x 356 mm (d) x 83 mm (h)
- Optional height rod available (item # 402KLROD)
- Optional wheels available (item # 55000)
- Optional counter weights available, to increase capacity to 490lb / 210kg (item # 55090)
- 10 Year Limited Warranty

**Optional Wheels**  
 Item # 55000



**Health o meter®**  
 P r o f e s s i o n a l  
[www.homscales.com](http://www.homscales.com)

Health o meter® Professional Scales  
 9500 West 55th Street USA  
 McCook, IL 60525-7110  
 Phone: 1-800-815-6615  
 Fax: 1-800-654-7330