

THE 499KL

WAIST-HIGH PRIMARY CARE DIGITAL SCALE

The 499KL is a value-engineered primary care scale. The scale provides the optimum in quality, reliability, and ease of use. Designed to meet the needs of doctor's offices, hospitals, clinics, acute and long-term care facilities, scale features a slim compact design, a large non-skid platform, and wheels that enable portability. The 499KL gives providers and patients a waist-high weight indicator for patients up to 460 pounds (209 kilograms), with a user-friendly digital interface that includes BMI, hold/release, zero, tare, and KG/LB conversion. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale

With its low platform height, the 499KL is easy to set up and offers patients a more comfortable and stable weighing experience.

499HB

"live" wrap-around handlebar

The 499HB is a unique handlebar accessory for the 499KL and is easy to attach to any new scale using the pre-drilled mounting holes and can also be attached to any existing scale using the easy-to-use drill guides. It provides the stability that many patients need without affecting the accuracy of the weight reading.



THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 499KL is reliable and easy to use.

120+ LOCAL
FIELD SPECIALISTS

#1 TRUSTED
BRAND

8" WIDE WHEELBASE
BEST IN ITS CLASS

100 YEARS
OF SERVICE

Health o meter®
Professional

Health o meter® Professional 499KL

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

460 lb / 209 kg

Resolution

0.2 lb / 0.1 kg

Display

1" / 25 mm, LCD display

Batteries

6 AA batteries included,
optional adapter
(order# ADPT31)

Functions

LB/KG conversion
LB/KG lock out
BMI
Zero
Tare
Hold/Release
Auto Zero
Auto Off
Parent/Child weighing function
Everlock® - permanent unit lock

Platform Size (w x d x h)

13 3/4" x 16 1/2" x 2 3/8"
(349 mm x 419 mm x 60 mm)
Dent and rust proof platform

Product Footprint

14 1/8" x 21 1/8" x 43 5/8"
(359 mm x 537 mm x 1108 mm)

Product Weight

21 lb / 9 kg

Portability

2 wheels, wide 8" wheelbase for ease
of movement

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
499KL	Waist-high digital scale, LB & KG
499KG	Waist-high digital scale, KG only
499KLHR	499KL with height rod included
499KGHR	499KG with height rod included
499KLAD	499KL with adapter ADPT31 included
499KGAD	499KG with adapter ADPT31 included
499KLROD	Optional height rod for 499 series
499HB	Optional handlebar accessory for 499 series scales
SS-499KL	ScaleSurance® 2 year extended warranty for 499 series*
ADPT31	Optional adapter
PLUGSET30	Plug adapters for use outside North America

499KLROD

Optional height rod



To learn more, call **1-800-815-6615**
or visit [homescapes.com/499KL](https://www.homescapes.com/499KL)

Health o meter®
Professional



*Not available in all countries.