

THE 524KL

DIGITAL PEDIATRIC SEAT SCALE

The 524KL pediatric scale was specifically created for nurseries, delivery rooms, and pediatric practices, making it an ideal choice for newborns, infants, and toddlers. With a maximum capacity of 50 pounds, this scale offers providers and patients a range of benefits, including safety, comfort, and flexibility. The deep infant seat on the 524KL ensures that the weight is easily and accurately captured while the seat supports the patient's head and body.

The 524KL display is positioned 17" above the scale base, ensuring an unobstructed easy-to-read weight measurement. The user-friendly digital interface includes features such as Zero, Tare, Pre-Tare, Hold/Release, and Auto-off, making it convenient and easy to use. Its motion detection technology ensures accurate weight measurements even when the patient moves.

The 524KL provides easy integration and seamless communication with today's EMR applications with the option to incorporate wireless connectivity for interfacing with computers, EMR software, or other electronic devices.



THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 524KL is your choice for safety, comfort, and reliability.

120+ LOCAL
FIELD SPECIALISTS

30-60 SEC
AUTO WEIGHT HOLD

#1 TRUSTED
BRAND

100 YEARS
OF SERVICE

Health o meter®
Professional

Health o meter® Professional 524KL

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

50 lb / 23 kg

Resolution

0 - 20 lb / 0.2 oz; 20 - 50 lb /
0.5 oz; 0 - 9 kg / 5g; 9 - 23 kg / 10 g

EMR Connectivity

Via optional Pelstar®
wireless technology

Display

1" / 25 mm, LCD display

Power Source

6 AA batteries included,
100-240V adapter (ADPT31)
not included

Functions

LB / KG conversion
LB / KG lock out
Zero
Tare
Pre-Tare
Hold / release
Reweigh
Recall
Auto zero
Auto off
Motion detection
Audible/mute sound option
Variable auto-off time
Enabling unit view with unit of measure locked
Auto-weight hold for 30 to 60 seconds
Everlock® - permanent weight lock

Seat Size (w x d x h)

24 ¼" x 14 ¼" x 18 ½"
(616 mm x 362 mm x 470 mm)

Product Footprint

24 ½" x 21" x 23 ⅛"
(622 mm x 533 mm x 587 mm)

Product Weight

21 lb / 10 kg

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
524KL	Digital Pediatric Seat Scale, LB & KG
SS-524KL	ScaleSurance® 2 year extended warranty for 524 series*
ADPT31	Optional adapter for 524 series of scales
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit homscales.com/ConnectivitySolutions



To learn more, call **1-800-815-6615**
or visit homscales.com/524KL

Health o meter®
Professional



*Not available in all countries.