THE 553KL

DIGITAL PEDIATRIC SCALE WITH EXTRA-WIDE TRAY

The 553KL pediatric scale was specifically created for nurseries, delivery rooms, and pediatric practices, making it an ideal choice for newborns, infants, and toddlers. With a maximum capacity of 44 pounds and 20 kilograms, this scale offers providers and patients a range of benefits, including safety, comfort, and flexibility. The scale ensures a secure

platform, keeps patients calm for accurate readings, and accommodates infants lying down to toddlers sitting up.

The 553KL pediatric scale offers a built-in measuring tape, allowing users to measure weight and height simultaneously. It also provides the flexibility of choosing between an optional digital or mechanical height rod. For easy portability and storage, there is an optional carrying case.

With its user-friendly digital interface, the 553KL scale includes features such as Zero, Tare, Pre-Tare, Hold/Release, Reweigh, and Auto-off, making it convenient and easy to use. Its motion detection technology ensures accurate weight measurements even when the patient moves.

The 553KL provides easy integration and seamless communication with today's EMR applications with the option to incorporate wireless connectivity for interfacing with computers, EMR software, or other electronic devices.

553KL-EHR

Scale with attachable digital height rod

THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 553KL is your choice for safety, comfort, and reliability.

120+ LOCAL FIELD SPECIALISTS

30-60 SEC AUTO WEIGHT HOLD

#1 TRUSTED
BRAND

100 YEARS OF SERVICE

weigheasier[®]



Health o meter® Professional 553KL

TECHNICAL SPECIFICATIONS

Capacity

44 lb / 20 kg

Resolution

0.5 oz / 10g

EMR Connectivity

Via optional Pelstar® wireless technology

Display

1"/ 25 mm, LCD display

Power Source

6 AA batteries included, 100-240V adapter (ADPT31)

not included

Functions

LB / KG conversion

LB / KG lock out

Zero Tare

Pre-Tare

Hold / release

Reweigh Recall

Auto zero Auto off

Motion detection

Variable auto-off time

Enabling unit view with unit of measure locked Auto-weight hold for 30 to 60 seconds

Everlock® - permanent weight lock

Tray Size ($w \times d \times h$)

24 %" x 13" x 1 %"

(632 mm x 330 mm x 35mm)

Product Footprint

24 %" x 16 ½" x 3 %"

(632 mm x 419 mm x 98 mm)

Product Weight

10 lb / 5 kg

Measuring Tape

0" - 24" / 0 cm - 58 cm

Portability

553CASE - optional carrying case

2210CART- optional rolling cart

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
553KL	Digital Pediatric Scale with Extra-Wide Tray, LB & KG
553KG	Digital Pediatric Scale with Extra-Wide Tray, KG only
553KL-HR	553KL with attachable mechanical height rod
553KL-EHR	553KL with attachable digital height rod
553HR	Mechanical height rod for 553 series of scales
553EHR	Digital height rod for 553 series of scales
SS-553KL	ScaleSurance® 2 year extended warranty for 553 series*
2210CART	Rolling scale/ medical device cart for 553 series of scales
553CASE	Carrying case for 553 series of scales
ADPT31	Optional adapter for 553 series of scales
ADPT30WA	Optional adapter and plugset for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit homscales.com/ConnectivitySolutions



2210CARTMedical Device Cart for 553
Series of Scales



Carrying Case for 553 Series of Scales

553CASE

To learn more, call 1-800-815-6615 or visit homscales.com/553KL



