


Stainless Steel* IV Pole with 4-Hook Top



MODEL: IVS-314

- Stainless steel* pole and base
- 4-hook detachable top
- Knob lock height adjustment
- Heavy duty, stainless steel base with low center of gravity design
- Ram's horn hooks
- 2" rubber wheel, ball bearing casters
- See available options on next page
- 45 lbs. load capacity under normal use

*Clinton stainless steel items are not M.R.I. compatible.

Shipping Weight: 17 lbs.

SPECIFICATIONS

Base Diameter—22"

Height Range: 60"—100"

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

OPTIONS:



IV-46 Pump Tray Support



IV-55 Extra Hooks and Organizer Hooks



IV-41 Universal Clamp



IV-50 Pole Outlet Strip



Heavy Duty Black Epoxy Wire Basket

IV-51B – 14" x 8" x 8.5"

IV-52B – 12" x 6" x 6.5"



Heavy Duty Stainless Steel Wire Basket

IV-51S – 14" x 8" x 8.5"

IV-52S – 12" x 6" x 6.5"



IV-57 Single Ram's Horn



IV-43 Steering Handle



IV-56 Drainage Hook



IV-44 Oxygen Cylinder Holder