

## Stainless Steel\*, Six-Leg, 4-Hook IV Pole



### MODEL: IVS-7334

- 4-hook top
- Large diameter, 6-leg epoxy-coated base
- Detachable ram's horn hooks with plastic end tips
- 3" rubber wheel ball bearing casters
- Stainless steel\* pole/knob-lock height adjustment
- See available options on next page
- 45 lbs. load capacity under normal use

\*Clinton stainless steel items are not M.R.I. compatible.

*Shipping Weight: 28 lbs.*

### SPECIFICATIONS

Base Diameter—26"

Height Range: 64"—100"

California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

## OPTIONS:



**IV-46 Pump Tray Support**



**IV-41 Universal Clamp**



**IV-42 Infusion Pump Frame**



**IV-50 Pole Outlet Strip**



**Heavy Duty Black Epoxy Wire Basket**

**IV-51B – 14" x 8" x 8.5"**

**IV-52B – 12" x 6" x 6.5"**



**Heavy Duty Stainless Steel Wire Basket**

**IV-51S – 14" x 8" x 8.5"**

**IV-52S – 12" x 6" x 6.5"**




**IV-57 Single Ram's Horn**



**IV-44 Oxygen Cylinder Holder**



**IV-56 Drainage Hook**



**IV-43 Steering Handle**



**IV-55 Extra Hooks and Organizer Hooks**