

MAICO Sessions

Let Your Data Speak Clearly

Fast and intuitive data management and visualization for audiometry & tympanometry tests

Master your Data, Focus on Care

Sessions is designed for fast and intuitive operation, allowing you to focus your attention on the patient. The simple navigation and modern design lets you grasp key details at a glance.

Simple Data Transfer

Your data is displayed the same way it is on the device, ensuring a consistent visual appearance.

Results transfer quick and easy with one click or automatically depending on the device.

Designed for Screening

The clean icon-based user interface enables an efficient workflow and high-speed data management. You can easily customize and print reports to guarantee a professional appearance of your business.

Fit for the Future

Sessions is a long-term solution serving as hub for future MAICO products. Your organization will benefit from MAICO Sessions for years to come.

Comprehensive Connectivity

You can link up Sessions to OtoAccess® or Noah databases with full functionality or install a standalone version for printing. Whatever your needs for data handling, Sessions provides a strong platform to work on.

Features at a Glance

- One click export as PDF or XML
- Standalone or database integration
- Simple and modern user interface
- Icon-based, intuitive navigation
- No training required
- Customizable reports



Technical Data Sessions

SUPPORTED DEVICES

MA 25e
MA 27e
MA 28
easyTymp
touchTymp and RaceCar

SUPPORTED MEASUREMENTS

Audiometry
Tympanometry and Acoustic Reflexes
Eustachian Tube Function (ETF) Intact and Perforated
Reflex Decay

TYPES OF EXPORT

PDF
XML

SUPPORTED DATABASES

OtoAccess®
Noah 4

PC-REQUIREMENTS

2 GHz Intel Core 2 Duo CPU
2 GB RAM
1 GB available disk space
USB drive

DISPLAY REQUIREMENTS

1024 x 768 resolution
Hardware accelerated DirectX/Direct3D graphics card


SYSTEM REQUIREMENTS

Windows® 10 SP1 (x86 and x64)
Windows® 8 / 8.1 (x86 and x64)
Windows® 7 SP1 (x86 and x64)
Preinstalled .NET Framework 4.5.2 or later

LANGUAGES

English, German, Spanish, French, Polish, Chinese


Specifications are subject to change without notice.



MAICO Sessions
PC software



OtoAccess®
Database



Noah Database