



## Hydraulic Transport Stretcher

Model NV8002



### Features:

- Fold-down side rails
- Six available I.V. receptacles
- Stationary foot section
- Six brake/steer pedals
- Dual pedestal hydraulics
- Dual side-mounted foot control
- Hands-free trendelenburg and reverse trendelenburg
- 4" Mattress included

### Specifications:

- Dimensions: 83 <sup>5</sup>/<sub>16</sub>" L X 33 <sup>9</sup>/<sub>16</sub>" W
- Weight Capacity: 750 lbs.
- Patient Surface: 75 <sup>7</sup>/<sub>16</sub>" L x 29 <sup>1</sup>/<sub>8</sub>" W
- Height Range: 22" - 35"
- Caster Diameter: 8"
- Positioning:
  - Head: 90°
  - Trendelenburg: 18°
  - Reverse Trendelenburg: 18°
- Warranty: 1 Year Parts & Labor

