



Stretchers

Model NV9000-OBGYN



Features:

- This OBGYN stretcher designed with functionality and convenience in mind, allowing caregivers to provide a more meaningful level of care to patients.
- One hand operated lowering leg section
- Deep well storage area for patient belongings
- Ergonomic design
- Built on a low height stretcher platform.
- Can serve as an interim transport stretcher for moving patients between departments.

Specifications:

- Two IV pole holders
- Roller bumpers
- Adjustable fowler section
- Fifth wheel steering
- Fold-down siderails
- 8" (20 cm) Tente® wheel casters
- Two side mounted pedals for patient surface height adjustment
- Four-caster braking system
- Ergonomic push handles
- Hydraulic lift system capable of supporting up to 550 lbs (250 kg)
- Head and foot-end mattress retainer
- Retractable water basin holder
- Overall Dimensions: 84" x 32.5"
- Surface Dimensions: 77.5" x 25" x 4"
- Trend/Reverse Trend: 15.5° - 17°
- Head Inclination: 81.5°
- Knee Inclination: 24°
- Height Range: 21.7" - 33"