

PowerSpin™ BX Centrifuge

Product Images



Short Description

POWERSPIN™ BX CENTRIFUGE
VARIABLE SPEED 1000-5000 RPM
8 PLACE FIXED ANGLE ROTOR
99 MIN DIGITAL TIMER
CAPACITY 8X10ML OR 8X12ML TUBES
POWER INPUT 110V

Description

UNICO is proud to introduce the latest addition to the PowerSpin™ Family – the BX model high speed centrifuge with 3 rotor options:

8-place Blood Tube rotor, 24 place Microhematocrit and 24 place Micro-Tube rotor

The state of the art design incorporates advanced safety features like Zero RPM Lid Lock and Auto Unbalance Detection, Rotor Recognition and Shutdown.

The built in oversized LCD Display provides real time data to the user. With the push of a button you can speed up or slow down the rotor within 1,000 - 11,000 RPM range (Max 5000 RPM with the Blood Tube Rotor). Maintenance free brushless DC motor is powerful and super quiet. The cover of the Microhematocrit rotor incorporates a reader, allowing for instant reading and eliminating the need for additional specimen handling.

Additional Information

Accessories	Header	Accessories
		C880-08 Fixed Angle Rotor, 8 place
		C880-02 Microhematocrit Rotor, 24 place
		C880-03 Microtube rotor, 24 place
		C880-10 BX Rotor Rack
		C881-03 Short tube adapter, for 13x75mm tubes
Specification		C881-04 Micro / Pediatric tube adapter for 0.5-2.0ml tubes
	Rotor	Fixed Angle, 8 place
	Speed Range	1000 -5000
	Maximum RCF(g)	2570
	Nominal Capacity	8 x 10ml
Size and Weight	Maximum Test Tube Size	17x105mm
	Size	10.5 in x 12.5 in x 9 in (270 mm x 320 mm x 220 mm)
	Shipping Size	14 in x 15 in x 12 in (350 mm x 375 mm x 300 mm)
	Weight	15.5 lbs (C881 & C882); 16.5 lbs (C883)
	Shipping Weight	19 lbs (C881 and C882) ; 20 lbs (C883)