

U-Base, Over Bed Table

MODEL: TS-165

- Gray laminate top
- Chrome plated base
- U-base for easy access
- Dual wheel nylon casters, two with locks
- Single hand, spring loaded height adjustments
- 50 lbs. load capacity (Distributed)

Shipping Weight: 29 lbs.

SPECIFICATIONS

Overall: 31.5" W x 16" D

Top: 31.5" W x 15" D

Base: 26.5" W x 16" D

Height Range: 31"—41.5"



California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1