


PowerSpin™ FX Centrifuge

Model: C808H

- POWERSPIN™ FX CENTRIFUGE
- FIXED SPEED 3400 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 6X10ML CAPACITY
- WITH 18 PLACE TUBE HOLDSTER™ RACK
- POWER INPUT 110V



Short Description

- POWERSPIN™ FX CENTRIFUGE
- FIXED SPEED 3400 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 6X10ML CAPACITY
- WITH 18 PLACE TUBE HOLDSTER™ RACK
- POWER INPUT 110V

UNICO® PowerSpin™ Series centrifuges are compact, economical and quiet. We offer both fixed (FX) and variable (LX) speed models with fixed angle rotors, for 6 or 8 tubes to spin 8– 15 ml tubes. By using adapters, you can also spin tubes ranging from 2 – 7 ml. An optional Tube Holdster™ allows you to have up to 18 tubes or adapters at your fingertips, side-mounted on the centrifuge. PowerSpins are pre-drilled/tapped to accept optional Tube Holdsters.

ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	Accessories
	C800-00 8 Place Rotor
	C800-01 8 Place Rotor
	C800-75 Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case
	C800-02 Regular tube shield, pack of 2
	C800-03 Adapter for short tubes, pack of 10
	C800-04 Pediatric or microtube adapter, pack of 2



	C800-05	Tube Lift Adapters, Pk/12
	C800-12	Balance Tube Holdster
SPECIFICATION	Number of Tubes	6
	Tube Capacity	2-15 ml
	Blood	Yes
	Specimens Fecal	Not recommended
	Urine	Not recommended
	Microhematocrit	Not recommended
	Speed	Fixed
	Rotor	Fixed Angle
	Maximum RPM	3400
	Maximum RCF(g)	1320
SIZE AND WEIGHT	Size	10.2 in x 11.8 in x 7.9 in
	Shipping Size	13.4 in x 14.8 in x 10.6 in
	Weight	13 lbs
	Shipping Weight	15 lbs

