


# PowerSpin™ PX Centrifuge

Model: C827

- POWERSPIN™ PORTA-SPIN™ CENTRIFUGE
- VARIABLE SPEED 500-3800 RPM
- 6 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 6X10ML CAPACITY
- POWER INPUT 15VDC



## Short Description

- POWERSPIN™ PORTA-SPIN™ CENTRIFUGE
- VARIABLE SPEED 500-3800 RPM
- 6 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 6X10ML CAPACITY
- POWER INPUT 15VDC

The C827 PXV 12V Variable Speed Centrifuge is a versatile and high-performance solution for mobile and field applications. Designed for 12V DC operation, it includes a vehicle power adapter (car adapter) for convenient use in remote locations, on-site testing, or mobile laboratories. This centrifuge offers variable speed control, allowing you to adjust the speed according to your specific sample processing needs. It features a 6-place rotor, accommodating up to six tubes for efficient sample separation. The dual lid lock ensures enhanced safety by securely locking the lid in place during operation, while the built-in leveling system provides stability even on uneven surfaces, ensuring consistent performance. With the C827 PXV, you can count on reliable and efficient centrifugation wherever you go, combining convenience and performance for a wide range of applications.



## ADDITIONAL INFORMATION

<b>TOP SELECTION</b>	No
<b>ACCESSORIES</b>	<b>Accessories (Sold Separately)</b> C800-00 8 Place Rotor



	C800-01	8 Place Rotor
	C800-75	Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case
	C800-02	Regular tube shield, pack of 2
	C800-03	Adapter for short tubes, pack of 10
	C800-04	Pediatric or microtube adapter, pack of 2
	C800-05	Tube Lift Adapters, Pk/12
	C800-12	Balance Tube Holdster
	C826-06	Porta-Spin 12 VDC Adapter, Plugs in to Car Charger, Replacement
	C826-12	Porta-Spin Ac Adapter, 110 Volt, USA Power Cord
<b>SPECIFICATION</b>	<b>Number of Tubes</b>	6
	<b>Tube Capacity</b>	2-15 ml
	<b>Blood</b>	Yes
	<b>Specimens Fecal</b>	Not recommended
	<b>Urine</b>	Not recommended
	<b>Microhematocrit</b>	Not recommended
	<b>Speed</b>	Variable
	<b>Rotor</b>	Fixed Angle
	<b>Maximum RPM</b>	3400
	<b>Maximum RCF(g)</b>	1320
<b>SIZE AND WEIGHT</b>	<b>Size</b>	10.2 in x 11.8 in x 7.9 in
	<b>Shipping Size</b>	13.4 in x 14.8 in x 10.6 in
	<b>Weight</b>	18 lbs
	<b>Shipping Weight</b>	19 lbs

