

PowerSpin™ PX Centrifuge

Model: C829

- POWERSPIN™ PORTA-SPIN™ CENTRIFUGE
- VARIABLE SPEED 500-3800 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 8X10ML CAPACITY
- POWER INPUT 15VDC



Short Description

- POWERSPIN™ PORTA-SPIN™ CENTRIFUGE
- VARIABLE SPEED 500-3800 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 8X10ML CAPACITY
- POWER INPUT 15VDC

The PortaSpin™ PX C829 Centrifuge is a robust and portable centrifuge designed to deliver reliable performance in mobile and field settings. Powered by a 12V DC vehicle power adapter (included), this centrifuge allows for easy operation directly from your vehicle's power source, making it perfect for remote locations, fieldwork, or mobile laboratories. With an 8-place rotor, the C829 PX accommodates up to eight tubes, enabling efficient processing of multiple samples at once. For added safety, it features a dual lid lock to ensure that the lid remains securely closed during operation. The built-in leveling system provides stability, ensuring consistent and reliable performance even on uneven surfaces. Whether you are conducting research in the field or need a centrifuge for on-the-go testing, the C828 PX offers the ideal solution, combining convenience, safety, and high performance in one portable unit.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	Accessories (Sold Separately) C800-00 8 Place Rotor



	C800-01	8 Place Rotor
	C800-75	Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case
	C800-02	Regular tube shield, pack of 2
	C800-03	Adapter for short tubes, pack of 10
	C800-04	Pediatric or microtube adapter, pack of 2
	C800-05	Tube Lift Adapters, Pk/12
	C800-12	Balance Tube Holdster
	C826-06	Porta-Spin 12 VDC Adapter, Plugs in to Car Charger, Replacement
	C826-12	Porta-Spin Ac Adapter, 110 Volt, USA Power Cord
SPECIFICATION	Number of Tubes	8
	Tube Capacity	2-15 ml
	Blood	Yes
	Specimens Fecal	Not recommended
	Urine	Not recommended
	Microhematocrit	Not recommended
	Speed	Variable
	Rotor	Fixed Angle
	Maximum RPM	3400
	Maximum RCF(g)	1320
SIZE AND WEIGHT	Size	10.2 in x 11.8 in x 7.9 in
	Shipping Size	13.4 in x 14.8 in x 10.6 in
	Weight	18 lbs
	Shipping Weight	19 lbs

