

Versa™ Fixed Height Exam Table

Specifications

Highly Customizable

- Easy-connect, removable upholstery top is available in standard & pediatric configurations
- Add options as needed to economically meet practice requirements

Improved Patient Experience

- Large, sturdy step with integral, anti-slip surface enhances patient safety
- Stable, cushioned platform with 500-lb patient weight capacity comfortably supports patients of all builds
- Optional multi-position patient-assist handles sized for more confident transitions and easier patient self-positioning
- Seamless upholstery & integral step surface offer easy cleanability and infection control

More Ergonomic Access

- Rounded corners and softened edges throughout safeguard patients and caregivers, while improving patient access
- Pass-through side drawer purposefully positioned at waist height to reduce bending and twisting
- Front drawer strategically positioned so caregivers do not infringe on a patient's personal space

Table Specifications

- **Weight Capacity:** 500 lbs (227 kg)
- **Seat height:** 32.5" (82.6 cm)
- **Length:** 58"; expandable to 72" (147.3 cm - 182.9 cm)
- **Pediatric table length:** 51"; expandable to 65" (129.5 cm - 165.1 cm)
- **Leg extension length:** 14" (35.6 cm)
- **Upholstery width:** 27.25" (69.2 cm)
- **Pediatric table width:** 26" (66 cm)
- 4.2 cubic feet (119 liters) of storage
- **Front drawer:** 20.9" x 18.8" x 3.8" (53.1 cm x 47.8 cm x 9.7 cm)
- **Side drawer:** 23" x 34.8" x 7" (58.4 cm x 88.4 cm x 17.8 cm)
- Paper roll storage capacity: accommodates up to 21" paper rolls
- Easy connect, removable upholstery top
- Industry-best 5 year warranty
- Easy-to-clean, seamless upholstery
- Two side drawer dividers for organization
- Large, sturdy patient step with integral anti-slip surface

Optional Features

- Hospital-grade electric receptacle
- Paper fastener strap on standard top offers convenience
- Patient assist handles
- 120VAC, 60Hz outlet
- 60-Watt front drawer heater
- 5° pelvic tilt
- Ergonomic Stirrups
- Footrest pad

Certifications

All Brewer products conform to US & Canadian electrical & medical industry codes & regulations and meet or exceed industry requirements.

- UL 60601-1
- CAN-CSA C22.2 No. 601.1
- ISO 13485:2016



Models:

Recommended Models:

- **1000AX-XXDX-XXX-__ [2 digit color code]**
Standard Model
- **1000AX-XSDT-RHX-__ [2 digit color code]**
Standard model with stirrups, pelvic tilt, receptacle and front drawer heater
- **1000AX-HXDX-XXX-__ [2 digit color code]**
Standard model with assist handles
- **1000AX-HSDT-RHX-__ [2 digit color code]**
Standard model with assist handles, stirrups, pelvic tilt, receptacle and drawer heater
- **1000PX-XXDX-XXX-__ [2 digit color code]**
Pediatric model


Accessories:


- **98464:**
Stainless Steel Treatment Pan
- **3500286:**
Deluxe articulating knee crutches with assembly kit
- **100253-TG:**
Urology Drain Pan System
- **3500288:**
Front Drawer Dividers
- **3500347:**
Brewer Welch Allyn® Exam Light / IV Bracket (factory installed)
- **3500398:**
Brewer Welch Allyn® Exam Light / IV Bracket (field installed)
- **3500305:**
Paper Fastener Strap
- **3500077:**
Footrest Pad


Create an elevated patient-provider experience.


Recommended Models	1000AX-XXDX-XXX	1000AX-XSDT-RHX	1000AX-HXDX-XXX	1000AX-HSDT-RHX	1000PX-XXDX-XXX
Weight Capacity	500	500	500	500	500
Standard or Pediatric	Standard	Standard	Standard	Standard	Pediatric
Patient Assist Handles			☑	☑	
Stirrups		☑		☑	
Drawer Warmer		☑		☑	
Pelvic Tilt		☑		☑	
Storage Capacity	4.2 cu ft	4.2 cu ft	4.2 cu ft	4.2 cu ft	4.2 cu ft
Side Drawer Dividers	☑	☑	☑	☑	☑
Footrest Pad	Optional	Optional	Optional	Optional	Optional


Available Colors.


 Clamshell
21


 Gunmetal
22


 Feather
23


 Black Satin
24


 Saddle
25


 Tapestry
26


 Cabernet
27

 Cocoa
28

 Deep Sea
29

 Ivy
30

 Deep Sapphire
31

 Blue Fog
32

Need a closer look? Contact us to receive physical samples of our vinyl or base materials.
Custom colors available.



1-888-273-9371 | BrewerCompany.com

The Brewer Difference.

Industry's most safety features.

Industry's only 5-year warranty.

USA-made in the Midwest.