

**Welch Allyn  
MicroTymp 4 Portable  
Tympanometer**

## Verify otoscopy findings and empower your diagnosis of middle ear conditions with objective evidence.

The **Welch Allyn MicroTymp 4** helps you verify middle ear conditions and improve diagnostic accuracy by providing objective analysis of tympanic membrane mobility and condition of the middle ear.

### FAST, FOR ALL PATIENTS


Test each ear in as little as three seconds from six months of age to geriatric patients.

### EFFICIENT AND COMPREHENSIVE EXAMS

Four button navigation to test L, R or both ears with automatic test measurements.



## USE THE MICROTYPMP 4 WHEN ADDRESSING COMPLAINTS OF HEARING LOSS.



Included TympView software provides a point for patient and parent education and helps support referrals. Results can be viewed on a PC and a PDF copy of the test results can be saved and printed.

Features two testing options, which automatically measure the following variables:

### 1 226 Hz TONE TYMPANOMETRY

- Tympanic peak pressure
- Admittance at peak
- Ear canal volume
- Gradient

### 2 ACOUSTIC REFLEX MEASUREMENT

- Ipsilateral stapedial reflexes
- Test using four frequencies
- Configurable reflex settings
- Automation detection

In addition to otitis media, the **MicroTymp 4** helps detect these other pathologies:

- Perforated tympanic membrane
- Patent tympanostomy
- Ossicular disruption
- Otosclerosis
- Tympanosclerosis
- Cholesteatoma

## TYMPANOMETRY TESTING HELPS SUPPORT ANTIBIOTIC STEWARDSHIP

Watchful waiting involves observing children with uncomplicated Acute Otitis Media (AOM) for 48 to 72 hours to allow for spontaneous recovery without antibiotics. The appropriate use of watchful waiting for children with AOM could reduce costs and improve health outcomes.\*

\*Lieberthal AS, Carroll AE, Chonmaitree T, et al. Clinical Practice Guideline: The Diagnosis and Management of Acute Otitis Media. Pediatrics. 2013;131(3):e964–e999



## TYMPANOMETRY CAN BE REIMBURSABLE USING THE FOLLOWING CPT CODE:<sup>\*</sup> 92567, TYMPANOMETRY (IMPEDANCE TESTING)

<sup>\*</sup>The information is provided for convenience only and represents no statement, promise or guarantee by **Welch Allyn** concerning coverage or levels of reimbursement. Payment will vary by geographic locality. It is always the provider's responsibility to determine accurate coding, coverage, and claim information for services provided.



## CONFIGURATIONS

Welch Allyn MicroTym 4 Portable Tympanometer	93700: System, Complete	93700-NP: System, No Printer
93710: MT4 Charging Cradle	●	●
93715: MT4 Power Supply	●	●
93720: Disposable MT4 Ear Tip Starter Set†	●	●
93730: MT4 Probe Tip Cleaning Kit	●	●
93750: MT4 Four-Volume Test Cavity	●	●
39414: MT4/OAE USB Data Cable	●	●
93771: MT4 Printer Cable	●	●
39410: MPT-II Thermal Printer Kit	●	—
93790-2: MT4 Manuals and Software USB Thumb Drive	●	●

For a complete list of **MicroTym 4** Disposable Ear Tips, refer to the Parts and Accessories tab on <https://www.hillrom.com/en/products/microtym-4-portable-tympanometric-instrument/>

## OPTIONAL ACCESSORIES AND REPLACEMENT PARTS

**93740: MicroTym 4 Probe Tip and Gasket Kit**

**93760: MicroTym 4 Portable Tympanometer Carry Case**

**39407: MPT-II Printer Power Supply/Battery Charger**

**39412: MPT-II Printer Paper Roll, 1 each**

**39416 : MPT-II Replacement Battery, Li-Ion**



Connect with your Hillrom sales representative today to learn more about our hearing screening portfolio.

**hillrom.com**

4341 State Street Road, Skaneateles Falls, NY 13153

Baxter, Hillrom, Welch Allyn and MicroTym are trademarks of Baxter International Inc. or its subsidiaries.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

APR277901 Rev 1 03-AUG-2022 ENG – US