

Unwrapped FAST
273°F / 134°C

9 min

Wrapped Class B
273°F / 134°C

15 min

Wrapped Class B
250°F / 121°C

45 min




Enbio PRO Autoclave Technical Data

Capacity



x16



Load space
11.2×7.5×2.8 in
(28.5 cm x 18 cm x 7 cm)

Max. packages amount (circa) **16**
Max. recommended batch weight **3.5 lb (1.6 kg)**
Process chamber volume **5.3 l**



Power supply

Enbio works with any typical electrical outlet.




Data archiving

All process data are automatically written to a USB disk.



Water supply

Enbio uses distilled water from external tank.




Process reports

Take full control of your sterilization process with **my.enbio app**.

Dimensions

22.6 in (57.5 cm)



10.6 in (27 cm)



7.9 in (20 cm)

Pressure & temperature

Working pressure **30.5 psi / 2.1 bar**
Max. pressure **33.4 psi / 2.45 bar**
Max. process temperature **278.6°F / 137°C**

Electrical

Max. power supply **110-120V / 60Hz**
Nominal power **3.25 kW**
Max. power consumption **15 A**

Misc.

Ingress protection **IP20**
Data archiving **USB**
Water conductivity* **< 15 µS/cm**

* - demineralized water only

Enbio Group AG
Eichengasse 3, 4702 Oensingen
Switzerland

US Distribution Center: Enbio Corp.
1075 Highway 124, Building 300, Suite 300
Hoschton, GA 30548



enbio PRO